

JW Multi Calibrator (Lot: 401) Target Value

JW BioScience Co.

2025.04 (Rev.01)

항 목	측 정 법	단 위	Conc.
A L B	BCG법	g/dl	3.5
A L P	PNPP/AMP/IFCC	IU/L	295
A M Y II	EPS/IFCC	IU/L	190
C A	Arsenazo III	mg/dl	10.2
T-CHO	효소법	mg/dl	165
C P K	IFCC법	IU/L	400
C R E A	Jaffe법	mg/dl	3.9
G G T	Glupa-C IFCC	IU/L	105
G L U(I)	GOD/POD법	mg/dl	190
G L U(II)	Hexokinase법	mg/dl	196
A S T	UV without P5P(IFCC법)	IU/L	114
A L T	UV without P5P(IFCC법)	IU/L	104
HDL-C	Direct법	mg/dl	54
I P	Phosphomolybdate, UV	mg/dl	5.3
LDL-C	Direct법	mg/dl	77
LDH	LD(P-L) Tris	IU/L	580
M G	Xylidyl Blue법	mg/dl	2.95
T G	Enzymetic	mg/dl	147
T P	Biuret법	g/dl	5.3
U A	효소법	mg/dl	5.45
B U N	Urease GLDH법	mg/dl	48.5
T-BIL (I)	Oxidation	mg/dl	4.1
D-BIL (I)	Oxidation	mg/dl	2.2
Lipase(III)	Colorimetric Method	IU/L	100

* 본 시약은 HA-7180에서 측정된 장비로서 기종에 따라 조금씩 오차가 발생할 수도 있음을 알려드립니다.

